

INSTRUCTIONS FOR USE OF MEDICAL DEVICE

FERTIMACY FEMININE COMFORT GEL™

Isotonic Gynaecological Solution, pH 4.0

SINGLE USE

3 ml pre-filled applicator

Composition:

Each 3 ml contains:

Sodium Edetate, Glycerine, Hydroxyethylcellulose, Potassium chloride, Lactic acid and Purified water.

Preservative - Polyhexamethylene biguanide hydrochloride (PHMB) 0,2 %.

Intended purpose:

FERTIMACY FEMININE COMFORT GEL is a feminine comfort gel used to relieve feminine discomfort, itching and odour. The action is local and solely to obtain or maintain the normal physiological conditions in the vagina.

Warnings or precautions:

The product must be administered exclusively by vaginal route. Products for topical use may cause sensitivity and should you be sensitive to the product, consult your doctor/pharmacist for an alternative treatment. Under normal use conditions, no overdosage harm has been reported.

Do not use the product if you are allergic to any of the components.

Do not use the product after the expiry date indicated on the package.

Do not use the product if the packaging is not sealed.

Do not swallow.

How to use this device:

Single use.

1. Remove the applicator from the protective wrapping
2. Remove the lid from the applicator.
3. Insert the applicator as deep as possible into the vagina
4. Slowly and gently push the plunger down to gently release the FERTIMACY FEMININE COMFORT GEL into the vagina until completely empty.
5. Remove the applicator from the vagina.
6. Discard the applicator after use.

Presentation:

3 ml pre-filled applicator

Storage:

Store at or below 25 °C. Protect from light.

Keep out of reach of children.

Name and business address of the holder of the certificate of registration

DUCIT

154 Ruimte road, Wierdapark, Centurion, South Africa

+27 12 654 3029

www.ducit.co.za

Date of issue or latest revision of the instructions for use: 2019 JAN